

## Error Strings

Error Code	Error Strings
10	"Slave absent"
20	"X position error"
21	"Y position error"
22	"Z position error"
23	"Homing error"
24	"X < minimum X"
25	"X > maximum X"
26	"Y < minimum Y"
27	"Y > maximum Y"
28	"Z < minimum Z"
29	"Z > maximum Z"
30	"X encoder fault"
31	"Y encoder fault"
32	"Z sensor fault"
33	"Safety activated"
34	"Invalid Z speed"
39	"Stop pressed"
99	"Unlisted motor error"
0	"Unknown error"

## External GSIOC Commands

Command	Type	Private	Description
%	I		Returns the character string: "223Vx.yz", where x, y, and z identify the software version for the sample changer.
\$	I		Returns a "\$" and resets the sample changer to its power-up state.
@	I		Reads non-volatile memory (NV-RAM) at current address. Returns "AA=xxxx" where:  AA - Value is the address (0 - 39). Xxxx - Data at the address.
@AA[=xxxx]	B		Sets the value at NV-RAM address where:  AA-Value is the address (0 - 39). xxxx-(Optional) Data at the address.
~n	BI		Sets test mode. Indicate one of the following for n:  1 - Performs XYZ test. <b>Must be run without needle from the home position.</b> 9 - Resets NV-RAM and initializes to defaults.

9	I	<p>Reads contact input event FIFO. If the queue is empty, “ 000000” is returned. If the queue is not empty, returns “Xttttt” where:</p> <p style="text-align: center;">X - State of the four contact inputs: 1 for closed, 0 for open. See table below.</p> <p style="text-align: center;">tttt - Time in 10 ms units since the last buffered 9 command.</p> <table><tr><td><i>if X=then</i></td><td>A=</td><td>B=</td><td>C=</td><td>D=</td></tr><tr><td>@</td><td>0</td><td>0</td><td>0</td><td>0</td></tr><tr><td>A</td><td>1</td><td>0</td><td>0</td><td>0</td></tr><tr><td>B</td><td>0</td><td>1</td><td>0</td><td>0</td></tr><tr><td>C</td><td>1</td><td>1</td><td>0</td><td>0</td></tr><tr><td>D</td><td>0</td><td>0</td><td>1</td><td>0</td></tr><tr><td>E</td><td>1</td><td>0</td><td>1</td><td>0</td></tr><tr><td>F</td><td>0</td><td>1</td><td>1</td><td>0</td></tr><tr><td>G</td><td>1</td><td>1</td><td>1</td><td>0</td></tr><tr><td>H</td><td>0</td><td>0</td><td>0</td><td>1</td></tr><tr><td>I</td><td>1</td><td>0</td><td>0</td><td>1</td></tr><tr><td>J</td><td>0</td><td>1</td><td>0</td><td>1</td></tr><tr><td>K</td><td>1</td><td>1</td><td>0</td><td>1</td></tr><tr><td>L</td><td>0</td><td>0</td><td>1</td><td>1</td></tr><tr><td>M</td><td>1</td><td>0</td><td>1</td><td>1</td></tr><tr><td>O</td><td>0</td><td>1</td><td>1</td><td>1</td></tr><tr><td>P</td><td>1</td><td>1</td><td>1</td><td>1</td></tr></table>	<i>if X=then</i>	A=	B=	C=	D=	@	0	0	0	0	A	1	0	0	0	B	0	1	0	0	C	1	1	0	0	D	0	0	1	0	E	1	0	1	0	F	0	1	1	0	G	1	1	1	0	H	0	0	0	1	I	1	0	0	1	J	0	1	0	1	K	1	1	0	1	L	0	0	1	1	M	1	0	1	1	O	0	1	1	1	P	1	1	1	1
<i>if X=then</i>	A=	B=	C=	D=																																																																																			
@	0	0	0	0																																																																																			
A	1	0	0	0																																																																																			
B	0	1	0	0																																																																																			
C	1	1	0	0																																																																																			
D	0	0	1	0																																																																																			
E	1	0	1	0																																																																																			
F	0	1	1	0																																																																																			
G	1	1	1	0																																																																																			
H	0	0	0	1																																																																																			
I	1	0	0	1																																																																																			
J	0	1	0	1																																																																																			
K	1	1	0	1																																																																																			
L	0	0	1	1																																																																																			
M	1	0	1	1																																																																																			
O	0	1	1	1																																																																																			
P	1	1	1	1																																																																																			
9	B	Clears the contact event FIFO and resets the contact input event timer.																																																																																					
a	I	Returns X travel length in Kilometers.																																																																																					
b	I	Returns Y travel length in Kilometers.																																																																																					
c	I	Returns Z travel length in Kilometers.																																																																																					
e	I	Reads the current error number. Returns “n” which identifies the error number; see page 5-1 for listing of errors. If no error has occurred, returns 0.																																																																																					
e	B	Clears error number.																																																																																					
Exyz H	B	<p>Sets X, Y, and Z motor status:</p> <p style="text-align: center;">X - 0 for disable or 1 for enable X motor. Y - 0 for disable or 1 for enable Y motor. Z - 0 for disable or 1 for enable Z motor.</p> <p>For example, the following command disables the motors: E000.</p>																																																																																					
H	B	Moves needle to home position.																																																																																					

## Gilson GSIOC Documentation

<b>I</b>	I		<p>Reads status of input contacts. Returns “cccc” where:</p> <p>cccc - Status of input contacts A, B, C, D: 1 if the contact is closed (shorted), 0 if open.</p>
<b>J</b>	I		<p>Reads status of output contacts and +12 V external auxiliary power. Returns “ccccp” where:</p> <p>cccc - Status of output contacts 1, 2, 3, and 4: 1 if the output is connected, 0 if disconnected.</p> <p>P - Status of +12 V external auxiliary power: 1 if connected, 0 if disconnected.</p> <p>Auxiliary power is off when sampler changer is turned on.</p>
<b>Jcccc[p]</b>	B		<p>Sets status of output contacts and +12 V external auxiliary power.</p> <p>Cccc - Output contacts 1, 2, 3, and 4: 1 to connect, 0 to disconnect, X for no change.</p> <p>P - (Optional) Auxiliary power: 1 to connect, 0 to disconnect, X for no change.</p>
<b>jc[ttt]</b>	B		<p>Pulses an output contact:</p> <p>c - Number of the output contact, 1 - 4. ttt - Duration of the pulse in tenths of seconds; default is 1.</p>
<b>Lx</b>	B		<p>Sets liquid level sensing threshold frequency based on current frequency and data at NV-RAM frequency threshold field.</p> <p>x - H for high frequency setting or L for low frequency.</p>
<b>M</b>	I		<p>Reads X, Y, Z motor status. Returns “xyz”. For each motor status, you see “U” for unpowered, “P” for powered, “R” for running, or “E” for error.</p>
<b>n</b>	I		<p>Reads the actual frequency of liquid level detector oscillator. Returns “ffff” which is frequency in Hz.</p>
<b>N</b>	I		<p>Reads the liquid level detector output. Returns “Iffff” where:</p> <p>I - A for air or L for liquid. ffff - Current sensitivity threshold frequency in Hz.</p>
<b>Nffff</b>	B		<p>Sets the liquid level sensing threshold frequency (ffff) in Hz.</p>

---

**Gilson GSIOC Documentation**

<b>P</b>	I		Reads the XYZ position in tenths of millimeters, based on the internal encoder wheel.
<b>S</b>	I		Reads the command in the synchronization buffer. Returns “ ” if buffer is empty.
<b>Smm</b>	B		Sends a synchronized buffered command (mm) that is executed when the sample changer is quiescent. Sending a command can overwrite unexecuted, existing commands. If you send this command without indicating a parameter (mm), the buffer is cleared.
<b>vzzzz,sss</b>	B		For tracking liquid height, raises or lowers the Z height of the needle at the designated speed:
<b>V</b>	I		Reads diverting valve status. Returns one of the following:  0 - Valve status off; the port connected to the needle is the one facing the rear of the sample changer.  1 - Valve status on; the port connected to the needle is the one facing you.
<b>Q</b>	I		Reads the Z travel range. Returns “min - max” where: